

Pressure Care

Seating 568 - 575

Supports 576

Elbow & Heel Protectors 577 - 578

Pressure Care **Classic-Med Cushion™**

The Classic-Med Cushion has a multi-stretch, vapour permeable and waterproof PU cover. Impregnated with acti-fresh biocide to help prevent the build up of bacteria, ensuring it can be easily cleaned to aid disinfection protocol.

1 Classic-Med Cushion™


Features

- Contoured cushion surface provides a greater weight distribution and a more comfortable support for the patient
- Internal air channels in the foam's surface help to reduce the build up of heat and moisture on the skin
- The cushion cover can be wiped down or completely removed and laundered up to 95°C to aid in decontamination protocol

091537356 430×430×100mm (17×17×4")

091537364 510×510×100mm (20×20×4")

Maximum user weight



Specification

Max. Weight Limit	17 stone / 108 kgs
Risk Factor	Low to Medium Risk
Warranty	2 years
Cushion Weight	2 kgs
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5



Memo-Pad Cushion™

The Memo-Pad Cushion has a multi-stretch, vapour permeable and waterproof PU welded cover. Impregnated with acti-fresh biocide to help prevent the build up of bacteria, ensuring it can be easily cleaned to aid disinfection protocol.

The high frequency welded seams and a fully concealed zip reduces the possibility of fluid ingress and aids infection control.

2 Memo-Pad Cushion™

Features


- The multi-stretch PU welded cover is waterproof and vapour permeable with acti-fresh biocide to help prevent the build up of bacteria ensuring, it can be easily cleaned to aid disinfection protocols
- The cushion cover can be wiped down or completely removed and laundered up to 95°C to aid decontamination protocols
- The temperature sensitive visco-elastic foam moulds to the shape of the sacral area providing greater weight distribution and pressure reduction

091537372 430×430×50mm (17×17×2")

091537380 430×430×76mm (17×17×3")


091537406 510×510×76mm (20×20×3")

Maximum user weight



Specification

Max. Weight Limit	23¾ stone / 152kgs
Risk Factor	Medium to High Risk
Warranty	2 years
Cushion Weight	3 kgs
Safety Standards	BS EN 597-1 Cig, BS EN 597-2 Match, BS7177: 1996 Crib Source 5



Apollo 4 and 4+ Cushion

The Apollo 4 is a high specification cushion that is ideally suited for use in an Acute Hospital, Nursing Home or Community environment also offering the Apollo 4+ while out and about still giving the client the same pressure reduction and care through the new innovation of the battery operated pump system that will run for up to 24 hours with an average client weight of 13 stone (82.7kgs).

If the pressure drops below the required setting on the Apollo4+, the system will automatically activate the compressors to maintain the support pressure by setting the cells to static mode.

The Apollo 4+ can also be supplied with a Patient hand controller which simply connects by USB to offer the client a little bit extra.

The waterproof vapour permeable cover on both models can be easily wiped clean or laundered at 95°C.

091537497 Apollo 4
091537505 Apollo 4+

Maximum
user weight



SUITABLE FOR THOSE AT
HIGH RISK OF DEVELOPING
PRESSURE ULCERS **HIGH**

Specification	
Max. Weight Limit	18 stone / 114.5kgs
Risk Factor	High Risk
Warranty	2 years
Cushion Dimensions	45×45×10
Cushion Weight	4 kgs
Power Unit Dimensions	Apollo 4 30×19×11cm Apollo 4+ 30×19×11cm
Power Unit Weight	2 kgs
Safety Standards	IEC601/1 BS5724 Part 1
Voltage	AC230v 50Hz
Noise Level	NC27
Power Input	12w
Fuse	T1A

APOLLO 4+ CUSHION



APOLLO 4 CUSHION



1 Homecraft Inflatable Seat Cushion

An inexpensive inflatable cushion for improved comfort and air circulation when seated. Inflates easily by mouth. Size 455×455mm (18×18"). Weight 200g.

091090968

Maximum user weight

113 kg **17³/₄ st**

SUITABLE FOR THOSE AT **LOW RISK** OF DEVELOPING PRESSURE ULCERS



2 Homecraft Eco Foam Cushion

This general purpose cushion is designed for basic pressure relief for wheelchair and seated users. Foam cushion with two way PU coated nylon stretch cover. Size 430×430×80mm (17×17×3¹/₄"). Weight 400g.

091091156

Maximum user weight

95 kg **15 st**

SUITABLE FOR THOSE AT **LOW RISK** OF DEVELOPING PRESSURE ULCERS



3 Homecraft Visco Wave Cushion

Unique visco wave upper surface, combined with high density foam that distributes pressure more evenly. The visco wave cushion has been designed for use in a wheelchair and/or day care chairs, by patients in the low/medium risk category. Supplied with a removable cotton cover. Size 420×420×75mm (16¹/₂×16¹/₂×3"). Weight 2kg.

091173566

Maximum user weight

75 kg **11³/₄ st**

SUITABLE FOR THOSE AT **LOW TO MEDIUM RISK** OF DEVELOPING PRESSURE ULCERS



4 Homecraft Posture Wave Cushion

The high resilient foam provides a comfortable upper surface with a central insert for posture support. Provides bony prominences protection at the base of the spine with integral safeguard foam inserts. The cushion has a removable cotton cover with a zip opening. Designed for use in a wheelchair and/or day care chairs, by patients in the low/medium risk category. Size 420×420×75mm (16¹/₂×16¹/₂×3"). Weight 2kg.

091173541

Maximum user weight

75 kg **11³/₄ st**

SUITABLE FOR THOSE AT **LOW TO MEDIUM RISK** OF DEVELOPING PRESSURE ULCERS



5 Homecraft Coccyx Wave Cushion

Superior wave upper surface with a high density foam insert, which offers greater cushioning to the coccyx. The vertex cushion provides coccyx protection at the base of the spine with an integral safeguard foam insert. Designed for use in a wheelchair and/or day care chairs, by patients in the low/medium risk category. Size 420×420×75mm (16¹/₂×16¹/₂×3"). Weight 2kg.

091173558

Maximum user weight

75 kg **11³/₄ st**

SUITABLE FOR THOSE AT **LOW TO MEDIUM RISK** OF DEVELOPING PRESSURE ULCERS



RETAIL PACKAGED
These products are supplied in retail packaging.



1 Homecraft Gel Plus Cushion

Development of the original Gel cushion has led to this unique Gel Plus cushion that combines combustion modified high resilient foam with visco elastic heat sensitive foam and silicon gel, to provide the user with the ultimate in comfort. Designed for those wheelchair users or seated patients that are at high risk of developing pressure sores. The cushion cover offers a PU coated 2 way stretch fabric top/side panels and non-slip base with sewn seams and zip opening. Size 430×430×80mm (17×17×3"). Weight 1.5kg
091163906

Maximum user weight **127 kg** **20 st** **SUITABLE FOR THOSE AT HIGH RISK OF DEVELOPING PRESSURE ULCERS** **HIGH**



1

2 Homecraft Visco Care Cushion

This Visco Elastic cushion offers affordable pressure relief for the high-risk user. Suitable for use in both wheelchairs and day care chairs, the anti-slip base provides additional security. The two-way stretch cover is waterproof and can be easily wiped clean. A water resistant zip allows the cover to be removed when required. Size 430×430 ×80mm (17×17×3¼"). Weight 1.1kg.
091091479

Maximum user weight **158 kg** **24¾ st** **SUITABLE FOR THOSE AT HIGH RISK OF DEVELOPING PRESSURE ULCERS** **HIGH**



2

3 Homecraft Gel Cushion

Wheelchair users and seated patients are at high risk of developing pressure sores. They need a reliable cushion that will prevent the onset of these sores and where necessary assist in the management of existing sores. This cushion has a carefully designed 'no lump' gel sac with bactericide incorporated into a foam exterior. It assists patient movement and gives excellent stability in use. It reduces interface pressures and provides high conformity giving specific ischial and perineal pressure protection. Covered in a breathable fabric that is waterproof and can be wiped clean. Two way PU coated nylon stretch permanent cover. Size 430×430×80mm (17×17× 3¼"). Weight 2.4kg.
091091164

Maximum user weight **108 kg** **17 st** **SUITABLE FOR THOSE AT HIGH RISK OF DEVELOPING PRESSURE ULCERS** **HIGH**



3

4 Homecraft Anyway Cushion

Many of the chairs in day rooms, wards and sitting areas are far from comfortable and certainly not pressure relieving. This general purpose cushion is designed to improve comfort and pressure distribution for wheelchair and seated users who have low to medium risk of pressure sores. It reduces interface pressures and can be used any way up and any way around. The multi-stretch PU coated nylon cover can be wiped clean easily. Size 480×480×80mm (19×19×3¼"). Weight 1kg.
091091172

Maximum user weight **98 kg** **15¼ st** **SUITABLE FOR THOSE AT LOW TO MEDIUM RISK OF DEVELOPING PRESSURE ULCERS** **LOW MID**




4

5 Homecraft Eco Plus Foam Cushion

Ideal for the medium risk user or larger person requiring a comfortable cushion when seated on a wheelchair or chair for a long period of time. This unique foam cushion can take up to a maximum user weight of 159kg (25st), with a height of 100mm (4") to provide that extra support. Two way stretch PU coated nylon permanent cover. Size 430×430×100mm (17×17×4"). Weight 1kg
091310754

Maximum user weight **159 kg** **25 st** **SUITABLE FOR THOSE AT MEDIUM RISK OF DEVELOPING PRESSURE ULCERS** **MID**



5



1 Dunlopillo Ring Cushion

Flame retardant bacteriostatic latex foam ring to give comfort when sitting. Can be steam sterilised. External diameter 432mm (17"). Internal diameter 127mm (5"). The cover is purchased separately. It is made from flame retardant cotton. Size 475×500mm (18³/₄×19³/₄").

Weight

091091008 **A** Cushion

1kg

091091032 **B** Dunlopillo Ring Cover

2 Visco Ring Cushion

An innovative visco foam ring that allows the user to sit for longer than is normally recommended for a ring cushion. The highly therapeutic material softens with the weight and heat of the body and moulds to the individual's requirements. The closed cell construction offers improved pressure distribution promoting increased sitting tolerance. An optional slip-over towelling cover is available. Size 430×430×100mm (17×17×4").

Weight

091091016 **A** Cushion 890g

Maximum user weight



3 Foam Ring Cushion

This durable foam cushion has a recess in the centre for comfort when seated. Extremely cost effective, it is supplied with a cover that can be removed for washing. Size 420×320mm (16¹/₂×12¹/₂"). Height 75mm (3"). Weight 350g.

091090869

4 Economy Rubber Ring Cushion

An inflatable rubber ring cushion to help provide additional comfort to the user. This ring cushion inflates easily by mouth, making it ideal to use when travelling. The surface of the ring is lightly patterned to prevent slipping. Use needle valve adaptor and a bicycle style pump to inflate. External diameter 457mm (18"). Internal diameter 160mm (6¹/₂"). Weight 715g.

091090984



1 Homecraft Inflatable PVC Ring Cushion

This ring cushion inflates easily by mouth, making it ideal for use when travelling. The surface has a suede like feel for additional comfort. External diameter 457mm (18"). Weight 150g. **091090992**

2 Commode Cushion

Clinically proven commode cushion that provides better pressure relief and distribution than visco-elastic and some leading air cushions. Washable at 40°C, the cushion can be used as supplied or with an optional autoclavable cover as illustrated. Size 500×500×100mm (20×20×4"). Weight 600g.

- 091081330** Cushion (White)
- 091081348** Cover Only (Turquoise)

Maximum user weight **90 kg** **14 st**

3 Commode Ring

This ring-shaped cushion can be used on most commodes or toilet seats to provide comfort. Covered in waterproofed cotton with a soft, fibre filling. External diameter 406mm (16"). Internal diameter 270mm (10½"). Weight 425g. **091081355**



homecraft



£ PRICES
For up-to-date prices please see our Website, alternatively contact your Sales Representative or Customer Services.



1 Polycore Wheelchair Seat Pad with Sides and Back

This seat pad with sides and back has been made to the Department of Health requirements and provides a high degree of support to patients suffering long-term immobility. Fully washable.

091438522

2 Heelift® Suspension Boot

The Heelift Suspension Boots relieve heel pressure by suspending the heel in space. It is ideal for pressure sore prevention and treatment, and for those with diabetes, cardiovascular disease, dietary deficiencies and more. Adjustable hook and loop straps and elevation pads allow for a custom fit; Tricot backing and polyethylene stiffener reduce friction and enhance client mobility. Prevents foot drop, unlike pillows or flat foam pads. The smooth Heelift is specifically designed for edematous legs to minimise chance for irritation or pressure points. Machine washable. Latex free.

Specifications

Smooth	Size	Calf Circumference	Weight Range
081422690	Petite	152-254mm (6-10")	32-59kg
081194398	Standard	203-356mm (8-14")	54-114kg
081422708	Bariatric	305-584mm (12-23")	100-272kg
Convuluted	Size	Calf Circumference	Weight Range
081433044	Petite	152-254mm (6-10")	32-59kg
081194380	Standard	203-356mm (8-14")	54-114kg



CONVOLUTED



SMOOTH

1 Wheelchair Fleece

This wool pile fleece brings warmth and comfort. It is machine washable, fire retardant and is secured with strong fabric ties to prevent slippage on the wheelchair. Size 890×432mm (35×17"). Weight 420g.

091091602

2 Padded Wheelchair Cushion

This quilted, fibre filled cushion, covers seat, back and arms of either chair or wheelchair. Seat size 360×260mm (14×10"). Weight 2kg.

091091487

3 Homecraft Curved Wheelchair Cushion

Firm crescent shaped foam pad to take sag out of a folding wheelchair seat. Size 432×432×51mm (17×17×2"). Weight 520g.

091091180

4 Homecraft Wheelchair Arm Covers

These simple but effective wheelchair arm covers fit snugly to the chair and give the wheelchair user more comfort. They are easy to attach with hook and loop fastenings and take up little space within the chair. Size 460×70×100mm (18×2¾×4"). Supplied as a pair. Weight 400g per pair.

091091198

5 Foam Seat Cushion

This high-grade foam cushion is ideal for a variety of situations. It can be used to elevate an easy chair or used in a wheelchair to provide additional comfort. Covered with a wipe clean PVC fabric. Size 430 x 430 x 80mm (17×17×3¼"). Weight 725g.

091090802


6 Luxury Fleece Topped Cushion

These Luxury Fleece Topped Cushions offer the utmost in comfort. They feature a luxury pure wool pile on the upper and front edge and have a vinyl underside. The cover is zipped so it can be removed for washing. Available in two sizes. Depth 406mm (16"). Height 50mm (2").

Width

091176338 406mm (16")

091176346 457mm (18")



1



1 Memory Foam Prop-up Wedge

Prop-up wedge manufactured with a memory foam comfort layer, offering optimum support whilst sitting or sleeping. Length 740mm (29¹/₂"). Width 580mm (23¹/₅"). Height 300mm (12"). Weight 4.8kg.

091164136

2 Memory Foam Seat Cushion

Memory foam seat cushion that can be used for extra support for existing chairs and is ideal when sitting for long periods. Length 440mm (17¹/₂"). Width 440mm (17¹/₂"). Height 50mm (2"). Weight 1.2kg.

091164128

3 Knee Support

This knee support made from polyurethane foam provides optimum comfort and pressure relief for the lower back and hips. It helps keep the back in the natural 'S' shape during sleep. Secured in place with a 'touch close' fastening. Covered with a soft polyester fleece. Available in visco elastic foam for optimum comfort, or high resilience foam for additional support. Dimensions 260×260mm (10¹/₄×10¹/₄").

091091636 High Resilience Foam

091091644 Visco Foam

4 Contour Leg Pillow

Provides comfort for the lower back, hips and knees. This leg pillow helps to align the legs, hips and spine into the correct position. Anatomically shaped to fit your legs comfortably all night. Supplied with a removable white cover. Length 250mm (10"). Width 200mm (8"). Depth between 70mm (2³/₄") and 150mm (5³/₄"). Weight 100g.

081258623

2



3



4



1 Rolyan® Elbow and Heel Protectors

Cushioned heel and elbow protectors with a foam pad which can be used in hospitals, nursing homes and extended care facilities. They can be machine washed and dried. They have a stretchable body and are knitted proportionately. Medication can be applied to pad when needed. One pair supplied.

	Size	Max Length	Weight
081288992	Small	230mm (9")	65g
081288984	Medium-Large	280mm (11")	70g
081289008	X-Large	330mm (13")	75g

2 Rolyan® Elbow and Heel Protectors with Shock Absorbing Padding

Helps protect the joint by diffusing and distributing pressure. Comfortable, stretchable stocking is proportionately knitted for a close fit. The shock absorbing Akton pads included with the protectors will not bottom out. May be used following a sudden injury or trauma to the elbow or heel, or simply to protect the joint as it rests against the surface. Pad is removable for washing. Sold as a pair.

	Size	Max Length	Weight
081288810	Small	230mm (9")	170g
081288828	Medium-Large	280mm (11")	175g
081288836	X-Large	330mm (13")	180g

3 Heel and Elbow Protectors

The 3mm thick polyurethane pad is held in place by a stretchable, soft terry loop material inside with a flat knit on the outside. The foam pad is tacked to the material to prevent slippage, yet allows the user to move the joint without a shearing effect. Five colour-coded adult sizes available. Machine washable. Sold in pairs.

	Size	Colour	Circumference Range	Length
081181718	Small	Yellow	150-240mm (6-9½")	235mm (9¼")
081181726	Medium	Blue	165-280mm (6½-11")	235mm (9¼")
081181734	Large	White	180-320mm (7-12½")	235mm (9¼")
081181742	X-Large	Green	190-343mm (7½-13½")	235mm (9¼")
081181759	XX-Large	Beige	230-420mm (9-16½")	310mm (12¼")

4 Deluxe Gel Heel and Elbow Protectors

Designed to improve comfort and help reduce shearing and friction on the heel or elbow. Features a removable gel pack and a plush terry loop material which is soft on the skin. Five colour-coded adult sizes available. Sold in pairs. To size, measure limb circumference.

	Size	Colour	Circumference Range	Length
081181767	Small	Yellow	150-240mm (6-9½")	235mm (9¼")
081181775	Medium	Blue	165-280mm (6½-11")	235mm (9¼")
081181783	Large	White	180-320mm (7-12½")	235mm (9¼")
081181791	X-Large	Green	190-343mm (7½-13½")	235mm (9¼")
081181809	XX-Large	Beige	230-420mm (9-16½")	310mm (12¼")





1 Fleece Heel Protectors

Washable wool pile pads that attach to the heel and fasten with padded straps. Provides comfort and minimises shearing tissue. Available to those with damaged or delicate heels. Supplied as a pair. Weight 210g.

091090455

2 Fleece Open Slippers

Comfortable slippers to protect the foot, ankle and heel. Made from lambswool with hook and loop straps. Supplied as a pair.

	Size	Weight
091090562	Small (Fits up to size 4)	160g
091090570	Medium (Fits up to size 5-9)	200g
091090588	Large (Fits up to size 10+)	270g

3 Gelbodies Heel and Elbow Protectors

The Trugel pressure relieving gel, enveloped in a soft and stretchy Coolband material, ensures these heel and elbow protectors give superior relief. Suitable for long term care, they can be used in the hospital and homecare environment. They are user friendly, lightweight, discreet and provide protection for both the heel and elbow from development of pressure sores. Sold in pairs. Latex free. 94% polyester, 6% spandex.

- Machine washable up to 40°C
- Gel pads can be cleaned with a damp cloth using a standard cleaning agent. The gel pads should not be machine washed
- Can be worn up to 12 hours in any 24 hour period

	Size	Dimensions
091090612	Small	200-250mm (8-10")
091090620	Medium	250-300mm (10-12")
091090638	Large	300-360mm (12-14")
091090646	X-Large	360-410mm (14-16")

4 Gelbodies Heel and Ankle Protectors

Made from the same high quality materials as the heel and elbow protectors above, but with a larger gel pad to protect the ankle bones and help prevent the development of pressure sores when side lying. Sold in pairs. Latex free. 94% polyester, 6% spandex.

- Machine washable up to 40°C
- Gel pads can be cleaned with a damp cloth using a standard cleaning agent. The gel pads should not be machine washed
- Can be worn up to 12 hours in any 24 hour period

	Size	Dimensions
091090653	Small	200-250mm (8-10")
091090661	Medium	250-300mm (10-12")
091090679	Large	300-360mm (12-14")

5 Gelbodies Universal Knee Protectors

Made from the same high quality materials as the Gelbodies protectors above, these universal knee protectors are designed to protect the knees from the development of pressure sores in the side lying position. Sold in pairs. Latex free. 94% polyester, 6% spandex.

- Machine washable up to 40°C
- Gel pads can be cleaned with a damp cloth using a standard cleaning agent. The gel pads should not be machine washed
- Can be worn up to 12 hours in any 24 hour period

	Size	Dimensions
091091511	Medium	280-330mm (11-13")
091091529	Large	330-410mm (13-16")

