

Falls Management

Fall Alarms	560 - 561
Patient Monitors & Alarms	562 - 563
Fall Mats	564
Home Safety	565
Days Red Walking Frames	566

Fall Alarms



How do I choose an Alarm?

Personal Magnet Alarms

Personal Magnet Alarms are easy to use, reliable and will suit just about every situation where a high fall risk patient needs monitoring. The pull cord clips to the patient's clothes and when the patient moves too far, the magnetic hold will disengage and sound the alarm.

Alarm Monitors

Alarm Monitors work in conjunction with sensor pads. They are good for low activity patients and only trigger the alarm if their body weight is removed from the sensor pad. Once the weight is back on the pad the alarm will go off. Can be used on beds or chairs and available in three versions, Economy, Standard and Deluxe.

ECONOMY

Low cost with basic functionality.

STANDARD

Tone switch, mains adaptor and nurse call jack.

DELUXE

All the same features as the standard with the added feature of voice record. Great for patients who simply need a reminder to sit or lay back down, call for help or who get agitated at the sound of an alarm beeping.



Magnet Alarm

This simple alarm is easy to operate, allows a quick response and is reliable and economical too. It is designed to be used on chairs with the cord attached to patient's clothes so that the alarm is activated when the cord is pulled as a result of the patient starting to stand. The pull string has a 25mm (1") diameter magnet on one end and a spring-loaded metal clip on the other. The alarm triggers and the red light illuminates when the magnet is pulled away from the alarm. A spring-loaded clip is included on the unit for easy attachment to a wheelchair.

Key features:

- The alarm unit measures 127×51×32mm (5×2×1¼"), making it small and discreet.
- Durable design – doesn't need an expensive rubber protector. Passed drop test by independent lab simulating 500 drops.
- A volume adjustment switch inside the battery compartment can be set at high (105dB) or low (90dB).
- The included 9V battery has an extra secure fit. When the alarm is dropped, the battery does not fall out of the alarm and also, makes it more difficult for patients to disable the alarm.
- Red, low battery light flashes when battery is low.
- Tool free battery door
- The length of the pull cord can be adjusted up to 762mm (30") allowing flexibility on how far the patient is able to move before activating the alarm.

- 081561489** Magnet Alarm
- 081561604** Replacement Magnet Pull Cord – Pack of 5

Alarm Monitors

For use with the chair sensor pads and bed sensor pads. The units are easy to move from wheelchairs to beds and the state-of-the-art technology ensures patient safety. These monitors, along with the pads, activate an alarm when a patient stands up or shifts weight off the pad.

- Sensor Pad Connect (green flash): Indicates that the alarm and pad have been connected and the alarm system is ready to use. The green flash goes off when the pad is lost, not connected or faulty.
- Low Battery Alert (amber flash): If the light is flashing, replace battery immediately.
- Alarm Alert (red flash): Patient has moved off the pad and needs assistance.
- Reset button needs to be pressed 3 times in 2 seconds to reset the alarm. This helps to prevent patient tampering.
- Includes a mounting strap, mounting hardware and a 9V battery.
- Economy and Standard Alarm Monitors have two tones to choose from and the Deluxe Alarm Monitor has one tone and a voice record feature.

- 081562537** Economy Alarm Monitor
- 081561497** Standard Alarm Monitor
- 081561505** Deluxe Alarm Monitor
- 091514520** Alarm Monitor Bed Bracket Clip
- 091527779** Mains Adaptor UK Plug



Code	Description	Connect Pad LED Indicator	Alarm LED Indicator	Low Battery LED Indicator	Reset Button	High/Low/Silent switch	0/1/2 second Delay Switch	Tone Switch	Voice Record	Nurse Call Jack	Mains Adaptor
081562537	Economy Alarm Monitor	x	x	x	x	x	x	x			
081561497	Standard Alarm Monitor	x	x	x	x	x	x	x		x	x
091514512	Deluxe Alarm Monitor	x	x	x	x	x	x	x	x	x	x

Alarm Sensor Pads

These durable vinyl-covered pads activate the alarm monitor when the patient rises or moves off the pad. Sensor pad set up instructions are printed on the pad for convenience along with a location to record the initial date of use.

Chair Sensor Pads

- 081561513** 254×381mm (10"×15") 90-Day
- 081561521** 254×381mm (10"×15") 180-Day
- 081561539** 254×381mm (10"×15") 1-Year

Bed Sensor Pads

- 081561547** 254×762mm (10"×30") 90-Day
- 081561562** 254×762mm (10"×30") 1-Year
- 081561588** 508×762mm (20"×30") 180-Day
- 081561596** 508×762mm (20"×30") 1-Year



Bed Sensor Pads



Chair Sensor Pads



90-Day Sensor Pads

Ideal for hospitals and patients who have a short-term stay

180-Day Sensor Pads

Useful for most rehab stays because of its durability

One-Year Sensor Pads

Long-term stays or at-home use benefit from this most durable pad

1



1 Fall Savers Floor Sentry Monitor

This versatile, easy-to-use, high quality bed and chair alarm system is an affordable tool for helping to manage fall risk patients and residents in care. The Fall Savers Monitor is a portable battery powered bed and chair alarm designed for use with either a bed or chair seat sensor pad. The monitor sounds an alarm when the patient's weight moves off the bed or chair seat sensor pad. Can also be used with a durable floor sensor pad to sound an alarm when the patient steps onto the floor sensor pad. The system is effective on patients weighing 28kg or more. With a high and low volume switch, the alarm can be adjusted to suit the current environment. To reduce false alarms, the monitor has an adjustable delay switch to assist the carer and patient. There is also a by-pass feature built-in, which allows the patient to be escorted out of a bed without sounding the alarm.

Floor pad 762×381mm (30½×5¼"). Compatible with SafePresence Sensor Mats.

091163054 Monitor

091162999 Floor Sensor Pad

2 Fall Savers Accessories

Not Illustrated

091163062 Nurse Call Connection Kit

3 SafePresence™ Bed Sensor Pad

The bed sensor is a thin flexible sensor pad used with the SafePresence Status Monitor and Audio Attendant, placed on top of the mattress and under the mattress cover without any discomfort to the person. The bed sensor detects the weight of the person being monitored and activates the monitor when the person leaves the bed. The robust sensor pad features antibacterial properties. Length 840mm (33½"). Width 500mm (20").

091310770

4 SafePresence™ Floor Sensor Pad

This heavy duty, floor sensor pad lies on the floor next to the bed or chair and connects to the SafePresence Audio Attendant monitor. When the person being monitored steps onto the floor sensor pad an alarm is activated on the Audio Attendant wireless remote. Length 890mm (35⅜"). Width 610mm (24⅔").

091310796

1



FLOOR SENSOR PAD

3



4









1 Helmet Skull Guard

Padded helmet that is flexible to comfortably conform to the contours of the head. Designed for children and adults with functional disorders that can make them liable to head injury, e.g. epilepsy, head drop syndrome, poor balance, autism, hydrocephalus, cerebral palsy with spastic diplegia, spinal cord hernia and haemophilia. The chin strap has a slide buckle to allow length adjustment.

To size, measure circumference of head.

	Size	Colour
091252907	46cm	● Bright Red
091252915	47cm	● Blue
091252923	48cm	● Blue
091252931	49cm	● Blue
091252949	50cm	● Blue
091252956	51cm	● Blue
091252964	52cm	● Blue
091252972	53cm	● Blue
091252980	54cm	● Blue
091252998	55cm	● Dark Brown
091253004	56cm	● Dark Brown
091253012	57cm	● Dark Brown
091253020	58cm	● Dark Brown
091253038	59cm	● Dark Brown
091253046	60cm	● Dark Brown
091253053	61cm	● Dark Brown



homecraft

2

homecraft

3

homecraft

4

homecraft

FOLDED

1 Homecraft Large Folding Bedside Mat

To be located alongside the bed for vulnerable patients that are at risk of falling. The mat will act as a cushion for the patient. Easily folded away, using the supplied securing strap and handle for transportation or storage. Can be easily wiped clean. Length 1700mm (67"). Width 890mm (35"). When folded length 850mm (33 1/2").

091173210

2 Homecraft Easy Access Bedside Mat

With a thickness of 50mm (2") and an angled edge, wheelchairs can easily roll on or off the mat. Made from high density foam with a vapour permeable cover with a non-slip base. The mat can be easily folded away, using the supplied securing strap and handle for transportation or storage. Can be easily wiped clean. Length 1700mm (67"). Width 650mm (25 1/2"). When folded length 850mm (33 1/2").

091173202

3 Homecraft Triple Folding Bedside Mat

A high density foam mat that is covered with a vapour permeable cover and a non-slip base. The mat can be easily folded away, using the supplied securing strap and handle for transportation or storage. Can be easily wiped clean. Length 1700mm (67"). Width 650mm (25 1/2"). When folded length 566mm (22").

091173194

4 Fall Out Mat

Suitable for use with a low level bed, this high density foam mat is covered in vinyl and folds for storage. Helps to reduce the possibility of impact related injuries when the situation dictates that bed rails should not be used. The mat folds into three for storage and transport. Size 1850x600x50mm (72 3/4x23 1/2x2").

091089630

1 Parnell Premier Cot Side Rail

The Parnell Premier Cot Side Rail has been designed for use with all domestic beds including adjustable beds. Tested by the University of Newcastle to BS EN60601-2-52:2010, the Parnell Premier Cot Side Rail helps to prevent patients from falling out of bed. Easy to assemble and fit with no tools, ties or straps required. The front footplates are flat so there is no tripping hazard or obstacle to wheelchairs or walking frames. Able to withstand a force/pull of 50kg in any horizontal direction when attached to the bed. Can be used with single, double or king size beds. Overall length 1220mm (48"). Overall height 870mm (34¹/₄"). The inner bars have a distance of 80mm (3¹/₄"). Minimum clearance required 30mm (1¹/₄"). Maximum clearance of 318mm (12¹/₂").

091322023

Maximum user weight

**2 Homecraft Cot Side Bumpers**

Made from thick foam, covered with vapour permeable, wipe clean fabric. These bumpers are made to fit a variety of cot sides and are easily fastened in place. Supplied as a pair. Length 1524mm (60"). Depth 381mm (15"). Weight 1kg each.

091089333**3 Dycem® Nose Over Toes Mat**

This innovative floor mat is designed to provide better traction underfoot and to prompt the user to position their feet correctly and adopt an appropriate position when transferring in or out of a chair.

091535541**4 Dycem® Non-Slip Floor Mat**

The Dycem Non-Slip Floor Mat has been designed to provide a stable, non-slip base for a patient when they are being transferred in any situation where the floor surface is slippery or potentially unsafe. Length 600mm (23¹/₂"). Width 450mm (18"). Weight 943g.





091102912 Blue**091102920** Grey**5 Home Safety Alert**

A wireless calling alarm system with one main receiver unit and two pendant calling transmitters. Call button and Panic button on each pendant. The main unit receives and recognises who is calling and the urgency using different sound and coloured light signals. A quick reset button on the main unit can cancel the call or alarm. Each transmitter uses 1 x 23A 12v battery (included) and the receiver unit is powered by a 6v adaptor (included) or 4 x AA batteries (not included).

091105105

RETAIL PACKAGED
These products are supplied in retail packaging.












Our most popular walking frames are now available in vibrant Metallic Red

These new RED Days Walking Frames were developed to meet demand from the NHS market for use with patients with dementia.

Bright colours can highlight important objects, making them easier for people with dementia to find and use making this range ideal to help increase the mobility for patients with dementia and reduce the risk of potential falls.

-  Wheeled and Non-Wheeled versions available
-  Hard anodised finish for durability
-  Comfortable latex free PVC handgrips
-  Low maintenance and easy to wipe clean
-  Lightweight aluminium frame
-  Simple push button height adjustment
-  Non-marking rubber ferrules

1 Days Adjustable Height Narrow Walking Frame

- Size**
- 091553890** Small
 - 091553866** Medium
 - 091553882** Large

2 Days Adjustable Height Standard Walking Frame

- Size**
- 091553775** Small
 - 091553767** Medium
 - 091553759** Large

3 Days Adjustable Height Narrow Wheeled Walking Frame

- Size**
- 091553874** Small
 - 091553841** Medium
 - 091553858** Large

4 Days Adjustable Height Standard Wheeled Walking Frame

- Size**
- 091553809** Small
 - 091553791** Medium
 - 091553783** Large

Spares

- 091356120** 25mm PVC ferrule
- 091553981** Wheeled Extension Legs

Specifications						
	Size	Overall Width B	Width at Front D	Depth C	Min. Height A	Max Height A
091553890	Small	59cm (23½")	41cm (16")	46cm (18")	70cm (27½")	78cm (30¾")
091553866	Medium	60cm (23¾")	42cm (16½")	46cm (18")	80.5cm (31¾")	88.5cm (34¾")
091553882	Large	62cm (24½")	42cm (16½")	46cm (18")	87.5 cm (34½")	95.5cm (37½")
091553775	Small	65cm (25½")	46cm (18")	49cm (19¼")	67cm (26½")	75cm (29½")
091553767	Medium	66cm (26")	47cm (18½")	55cm (21½")	77cm (30¼")	85cm (32½")
091553759	Large	67cm (26½")	51cm (20")	57cm (22½")	87cm (34¼")	95cm (37½")
091553874	Small	59cm (23½")	50cm (19¾")	49cm (19¼")	73cm (28¾")	81cm (31¾")
091553841	Medium	60cm (23¾")	51cm (20")	49cm (19¼")	83.5cm (33")	91.5cm (36")
091553858	Large	62cm (24½")	51cm (20")	49cm (19¼")	90.5cm (35¾")	98.5cm (38¾")
091553809	Small	65cm (25½")	55cm (21½")	52cm (20½")	70cm (27½")	78cm (30¾")
091553791	Medium	66cm (26")	56cm (22")	58cm (22¾")	81cm (31¾")	89cm (35")
091553783	Large	67cm (26½")	51cm (20")	57cm (22½")	87cm (34¼")	95cm (37½")